

# Evanesce™ flow

flowable restorative

from **CLINICIAN'S CHOICE**  
DENTAL PRODUCTS INC.

RESTORATIVE

COMPOSITE



## Today's minimally-invasive dentistry demands more from flowable composites.

In addition to its traditional use as a base liner under posterior composite restorations, the ability to provide a strong, durable, and esthetic conservative restoration using a flowable composite is an invaluable asset to any dental practice.

All-new **EVANESCE™ FLOW** has been engineered to complement the Evanesce Universal Restorative composite line and features the very same chemistry; resulting in amazing esthetics, high strength and preferred handling. The shades and opacities available in Evanesce Flow mirror the applications that Evanesce is designed for, and provide an excellent shade and opacity match to Evanesce and the VITA Classic shade guide, enabling it to be used on its own or in combination with Evanesce. Evanesce Flow moves easily when manipulated without sticking to the instrument causing pull-back. Evanesce Flow is thixotropic for improved control, is easily dispensed from the syringe, and will not slump on a vertical surface - making it ideal for use in Class V's and small Class II's. Regardless of its application, its relatively low rate of shrinkage (3.6% compared with up to 5%), and high filler content (62% by weight and 48% by volume), make it truly universal for use anywhere in the mouth.



**EVANESCE FLOW** is clearly distinguishable from composite on a radiograph for ease of diagnosis.

## Whether you're choosing for shade, location or opacity – there's an EVANESCE FLOW for that!

Evanesce Flow is ideal for superficial Class I and II's due to its esthetics, ability to be adapted to the tooth and high strength. Class V's benefit from its no-slump handling and range of opacities and shades. Evanesce Flow is invaluable as a liner base in posterior restorations to support the dentinal/adhesive seal, preventing post-op sensitivity due to its low contraction shrinkage. Radiographically, Evanesce Flow is easily distinguishable from tooth structure and neighbouring composite with a radiopacity of 200% Al.

Available in the most popular shades, Evanesce Flow A1, A2, and A3 Enamel shades will satisfy the majority of applications. However, an A4 Universal shade is also available for the aging population; while B1 Enamel and BL1 (Bleach) Universal shades are perfect for younger patients or for those who have previously whitened their teeth. In situations where there are varying degrees of cavity depth or discoloration of the underlying tooth structure, A4 Universal and A2 Dentin shades provide the ability to esthetically mask-out or layer the restoration.

Esthetic and strong, easy to express and manipulate, available in the most popular shades and opacities, Evanesce Flow is the ideal solution to minimally-invasive restorations and a vital component to your large composite cases.



### EVANESCE FLOW

1.7gm syringe

A1 Enamel	323501
A2 Enamel	323502
A3 Enamel	323503
B1 Enamel	323505
A2 Dentin	323511
A4 Universal	323520
B1 Universal	323525

### Single-dose 20pk

A1 Enamel	323531
A2 Enamel	323532
A3 Enamel	323533
B1 Enamel	323535
A2 Dentin	323546
A4 Universal	323550
B1 Universal	323555



TEACHING BETTER DENTISTRY

1-800-265-3444  
clinicalresearchdental.com